


Baxter

AK 98

DIALYSIS MACHINE



1. Venous pressure control
2. Arterial pressure control
3. Machine state indicator
4. Time indicator
5. Blood path
6. Fluid path
7. Blood pressure monitor (BPM) button
8. Blood pressure monitor (BPM) read out field
9. Diascan read out field
10. Diascan button
11. Treatment overview field
12. Treatment overview field
13. Treatment overview field
14. Treatment overview field
15. Treatment overview field
16. Alarm tab
17. Information tab
18. Treatment graph tab
19. Information field
20. Patient page
21. Functions button
22. Disinfection button
23. Blood button
24. Fluid button
25. Fluid bypass button
26. Ultrafiltration (UF) button
27. Treatment History page
28. Status bar
29. Monitor identifier (Nickname)
30. On/Off button
31. Schedule indicator
32. Main switch indicator
33. Mute button
34. Blood pump up button
35. Blood pump button
36. Blood pump down button

AK 98

Operator's Panel

Machine Introduction SW 3.XX

Important notice: The information in this document does not relieve any user from his or her duty to carefully read the full text of the Operator's Manual for the **AK 98** dialysis machine.

Baxter & AK 98 are trademarks of Baxter International Inc., or its subsidiaries.